

Rabbit Anti-Guinea Pig IgG Heavy & Light Chain antibody {Unconjugated} (STJ99572)

STJ99572

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Unconjugated rabbit polyclonal anti-Guinea Pig IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Rabbit
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	4-5 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
Formulation	PBS, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	