

Goat Anti-Mouse IgG Heavy & Light Chain antibody {Unconjugated} (STJ99671)

STJ99671

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Unconjugated goat polyclonal anti-Mouse IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB, IHC, IF, ICC, ELISA
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:500-1000 IF 1:5000-10000
Formulation	PBS, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	