

Anti-MHC 2 antibody [M5/114] {PE} (STJA0000054)

STJA0000054

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-MHC 2 (I-A/I-E) is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | M5/114 |
| Concentration | 0.2 mg/mL |
| Conjugation | PE |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2bk |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 14961 |
| Gene Symbol | H2-Ab1 |
| Uniprot ID | HB2A_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



Mouse splenocytes were stained with Anti-Mouse I-A/I-E Monoclonal Antibody (PE Conjugated) (Used at 0.02 µg/10⁶ cells dilution) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081