

Anti-MHC 2 antibody [M5/114] {AF488} (STJA0000060)

STJA0000060

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MHC 2 (I-A/I-E) is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5/114
Concentration	0.5 mg/mL
Conjugation	AF488
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14961
Gene Symbol	H2-Ab1
Uniprot ID	HB2A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse splenocytes were stained with Anti-Mouse I-A/I-E Monoclonal Antibody (AF488 Conjugated) (Used at 0.2 µg/10⁶ cells dilution) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081