

Anti-CD83 antibody [HB15e] {APC} (STJA0000093)

STJA0000093

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD83 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HB15e
Concentration	
Conjugation	APC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	9308
Gene Symbol	CD83
Uniprot ID	CD83_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



PMA and ionomycin-stimulated (4h) Jurkat cells were stained with Anti-CD83 antibody (STJA0000093) (filled gray histogram) or Mouse IgG1, κ Isotype Control (APC Conjugated) (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081