

Anti-CD47 antibody [MIAP410] {PerCP} (STJA0000202)

STJA0000202

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD47 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIAP410
Concentration	
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
Uniprot ID	CD47_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD47 antibody (STJA0000202) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.