

Anti-CD47 antibody [MIAP410] {PerCP/Cyanine5.5} (STJA0000212)

STJA0000212

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD47 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIAP410
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
Uniprot ID	CD47_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD47 antibody (STJA0000212) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081