

## Anti-CD62L/SELL antibody [DREG56] {FITC} (STJA0000344)

STJA0000344

### GENERAL INFORMATION


|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Mouse monoclonal antibody anti-L-selectin is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC   |
| <b>Host/Source</b>       | Mouse  |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | DREG56   |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | FITC   |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG1k  |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |             |
|---------------------------|-------------|
| <b>Gene ID</b>            | 6402        |
| <b>Gene Symbol</b>        | SELL        |
| <b>Uniprot ID</b>         | LYAM1_HUMAN |
| <b>Immunogen</b>          |             |
| <b>Immunogen Region</b>   |             |
| <b>Specificity</b>        |             |
| <b>Immunogen Sequence</b> |             |



Human peripheral blood lymphocytes were stained with Anti-CD62L antibody (STJA0000344) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081