

## Anti-CD147/BSG antibody [HIM6] {AF647} (STJA0000368)

STJA0000368

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Basigin is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC  |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | HIM6   |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | AF647  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG1k  |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            | 682        |
| <b>Gene Symbol</b>        | BSG        |
| <b>Uniprot ID</b>         | BASI_HUMAN |
| <b>Immunogen</b>          |            |
| <b>Immunogen Region</b>   |            |
| <b>Specificity</b>        |            |
| <b>Immunogen Sequence</b> |            |



Human peripheral blood lymphocytes were stained with Anti-CD147 antibody (STJA0000368) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081