

Anti-CD11b antibody [M1/70] {FITC} (STJA0000430)

STJA0000430

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/70
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Itgam
Uniprot ID	ITAM_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse bone marrow cells were stained with Anti-Mouse/Human CD11b antibody (STJA0000430) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081