

Anti-CD9 antibody [HI9a] {PE/Cyanine5} (STJA0000446)

STJA0000446

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD9 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HI9a
Concentration	
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	928
Gene Symbol	CD9
Uniprot ID	CD9_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Human platelets were stained with Anti-CD9 antibody (STJA0000446) (filled gray histogram). Unstained platelets (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081