

## Anti-CD90.2 antibody [30H12] {FITC} (STJA0000502)

STJA0000502

### GENERAL INFORMATION


|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rat monoclonal antibody anti-CD90.2 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC   |
| <b>Host/Source</b>       | Rat  |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | 30H12  |
| <b>Concentration</b>       | 0.5 mg/mL  |
| <b>Conjugation</b>         | FITC   |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG2bk   |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            | 21838      |
| <b>Gene Symbol</b>        | Thy1       |
| <b>Uniprot ID</b>         | THY1_MOUSE |
| <b>Immunogen</b>          |            |
| <b>Immunogen Region</b>   |            |
| <b>Specificity</b>        |            |
| <b>Immunogen Sequence</b> |            |



C57BL/6 murine splenocytes were stained with Anti-CD90.2 antibody (STJA0000502) [Used at 0.2 µg/10<sup>6</sup> cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081