

Anti-Cd4 antibody [GK1.5] {PerCP} (STJA0000532)

STJA0000532

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GK1.5
Concentration	0.2 mg/mL
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
Uniprot ID	CD4_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse splenocytes were stained with Anti-Mouse CD4 Monoclonal Antibody (PerCP Conjugated) (Used at 0.05 µg/10⁶ cells dilution) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081