

Anti-Cd8a antibody [OX-8] {PerCP/Cyanine5.5} (STJA0000542)

STJA0000542

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD8 alpha chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-8
Concentration	
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24930
Gene Symbol	Cd8a
Uniprot ID	CD8A_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Rat splenocytes were stained with Anti-CD8a antibody (STJA0000542) (filled gray histogram) or Mouse IgG1, κ isotype Control (PerCP/Cyanine5.5 Conjugated) (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081