

Anti-CD44 antibody [IM7] {AF488} (STJA0000554)

STJA0000554

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD44 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IM7
Concentration	
Conjugation	AF488
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12505
Gene Symbol	Cd44
Uniprot ID	CD44_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-Mouse/Human CD44 antibody (STJA0000554) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.