

Anti-CD25/Il2ra antibody [PC-61.5.3] {FITC} (STJA0000568)

STJA0000568

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PC-61.5.3
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
Uniprot ID	IL2RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD25 antibody (STJA0000568) [Used at 0.2 µg/10⁶ cells dilution] and Anti-Mouse CD4 Monoclonal Antibody (APC Conjugated) (right). Splenocytes stained with Anti-Mouse CD4 Monoclonal Antibody (APC Conjugated) are used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081