

Anti-Ly6g antibody [1A8] {FITC} (STJA0000623)

STJA0000623

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Lymphocyte antigen 6G is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 1A8 |
| Concentration | 0.5 mg/mL |
| Conjugation | FITC |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 546644 |
| Gene Symbol | Ly6g |
| Uniprot ID | LY6G_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



Mouse bone marrow cells were stained with Anti-Ly6G antibody (STJA0000623) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081