

Anti-Ly-6G/Ly-6C antibody [RB6-8C5] {PE/Cyanine7} (STJA0000731)

STJA0000731

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly-6G/Ly-6C (Gr-1) is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RB6-8C5
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	546644
Gene Symbol	Ly6g
Uniprot ID	LY6G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine bone marrow cells were stained with Anti-Mouse Gr-1 Monoclonal Antibody (PE/Cyanine7 Conjugated) (Used at 0.2 μ g/ 10^6 cells dilution) (filled gray histogram). Unstained bone marrow cells (empty gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081