

Anti-TER-119 antibody [TER-119] {PerCP/Cyanine5.5} (STJA0000795)

STJA0000795

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TER-119 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TER-119
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence



C57BL/6 murine bone marrow cells were stained with Anti-TER-119 antibody (STJA0000795) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.