

Anti-CD279/Pdcd1 antibody [29F.1A12] {PE/Cyanine5.5} (STJA0000831)

STJA0000831

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Programmed cell death protein 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29F.1A12
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	18566
Gene Symbol	Pdcd1
Uniprot ID	PDCD1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Con-A stimulated C57BL/6 splenocytes (3 days) were stained with Anti-CD279/PD-1 antibody (STJA0000831) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081