

Anti-Cd274/PD-L1 antibody [29E.2A3] {PE/Cyanine7} (STJA0000842)

STJA0000842

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Programmed cell death 1 ligand 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29E.2A3
Concentration	
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
Uniprot ID	PD1L1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



PHA-stimulated (3 days) human peripheral blood lymphocytes were stained with Anti-CD274/PD-L1 antibody (STJA0000842) (filled gray histogram) or Mouse IgG1 Isotype Control PE/Cy7 (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081