

Anti-CD69 antibody [FN50] {AF488} (STJA0000879)

STJA0000879

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Early activation antigen CD69 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FN50
Concentration	
Conjugation	AF488
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	969
Gene Symbol	CD69
Uniprot ID	CD69_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



PMA and ionomycin-stimulated (4h) Jurkat cells were stained with Anti-CD69 antibody (STJA0000879) (filled gray histogram) or Mouse IgG1, κ Isotype Control (AF488 Conjugated) (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081