

Anti-CD11b/ITGAM antibody [ICRF44] {FITC} (STJA0000909)

STJA0000909

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin alpha-M is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ICRF44
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3684
Gene Symbol	ITGAM
Uniprot ID	ITAM_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Human peripheral blood granulocytes were stained with Anti-CD11b antibody (STJA0000909) (filled gray histogram) or Mouse IgG1 Isotype Control FITC (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081