

Anti-CD107a/LAMP-1 antibody [H4A3] {PE} (STJA0000931)

STJA0000931

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lysosome-associated membrane glycoprotein 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H4A3
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3916
Gene Symbol	LAMP1
Uniprot ID	LAMP1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Intracellular staining of the Jurkat cell line with Anti-Human CD107a Monoclonal Antibody (PE Conjugated) (filled gray histogram) or Mouse IgG1 Isotype Control PE (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081