

Anti-CD90/THY1 antibody [5E10] {PerCP/Cyanine5.5} (STJA0000961)

STJA0000961

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Thy-1 membrane glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5E10
Concentration	
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	7070
Gene Symbol	THY1
Uniprot ID	THY1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Human erythroleukemic cell line (HEL) were stained with Anti-CD90 antibody (STJA0000961) (filled gray histogram) or Mouse IgG1, κ Isotype Control (PerCP/Cyanine5.5 Conjugated) (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081