

Anti-CD64/Fcgr1 antibody [X54-5/7.1] {APC} (STJA0001027)

STJA0001027

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-High affinity immunoglobulin gamma Fc receptor I is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	X54-5/7.1
Concentration	0.2 mg/mL
Conjugation	APC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14129
Gene Symbol	Fcgr1
Uniprot ID	FCGR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine bone marrow cells were stained with Anti-CD64 antibody (STJA0001027) [Used at 0.2 µg/10⁶ cells dilution] and Anti-Mouse CD11b Monoclonal Antibody (FITC Conjugated).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081