

Anti-CD44H antibody [OX-49] {PerCP/Cyanine5.5} (STJA0001073)

STJA0001073

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD44H is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-49
Concentration	
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence



Rat splenocytes were stained with Anti-CD44H antibody (STJA0001073) (filled gray histogram) or Mouse IgG2a, κ Isotype Control (PerCP/Cyanine5.5 Conjugated) (empty black histogram).