

Anti-CD90.1 antibody [OX-7] {FITC} (STJA0001082)

STJA0001082

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD90/CD90.1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-7
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24832
Gene Symbol	Thy1
Uniprot ID	THY1_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Rat thymocytes were stained with Anti-CD90/Mouse CD90.1 antibody (STJA0001082) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained thymocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081