

Anti-Cd19 antibody [1D3] (STJA0001211)

STJA0001211

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-lymphocyte antigen CD19 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D3
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12478
Gene Symbol	Cd19
Uniprot ID	CD19_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	