

Anti-MHC 2 (I-A/I-E) antibody [M5/114] (STJA0001219)

STJA0001219

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MHC 2 (I-A/I-E) is suitable for use in Flow Cytometry, Western Blot, Immunoprecipitation and Immunofluorescence research applications.
Applications	FC, WB, IP, IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5/114
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	14961
Gene Symbol	H2-Ab1
Uniprot ID	HB2A_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	