

Anti-MHC 2 (I-A/I-E) antibody [M5/114] {PE/Cyanine5.5} (STJA0001222)

STJA0001222

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MHC 2 (I-A/I-E) is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5/114
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14961
Gene Symbol	H2-Ab1
Uniprot ID	HB2A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	