

Anti-Cd86 antibody [GL-1] (STJA0001235)

STJA0001235

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-lymphocyte activation antigen CD86 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, IHC, IF, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GL-1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12524
Gene Symbol	Cd86
Uniprot ID	CD86_MOUSE
Immunogen	
Region	
Specificity	
Immunogen Sequence	