

Anti-CD16/32 antibody [2.4G2] (STJA0001242)

STJA0001242

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD16/32 is suitable for use in Universal and Flow Cytometry research applications.
Applications	Block, FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2.4G2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	14130
Gene Symbol	Fcgr2
Uniprot ID	FCGR2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	