

Anti-CD16/32 antibody [2.4G2] {PE/Cyanine7} (STJA0001246)

STJA0001246

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-CD16/32 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 2.4G2 |
| Concentration | |
| Conjugation | PE/Cyanine7 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2bk |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 14130 |
| Gene Symbol | Fcgr2 |
| Uniprot ID | FCGR2_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |