

Anti-Cd48 antibody [HM48-1] (STJA0001296)

STJA0001296

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-CD48 antigen is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, IHC, IP
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM48-1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12506
Gene Symbol	Cd48
Uniprot ID	CD48_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	