

Anti-CD54/Icam1 antibody [YN1/1.7.4] (STJA0001304)

STJA0001304

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Intercellular adhesion molecule 1 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC, IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	YN1/1.7.4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	15894
Gene Symbol	Icam1
Uniprot ID	ICAM1_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	