

Anti-CD54/Icam1 antibody [YN1/1.7.4] (STJA0001305)

STJA0001305

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Intercellular adhesion molecule 1 is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, WB, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	YN1/1.7.4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	15894
Gene Symbol	Icam1
Uniprot ID	ICAM1_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	