

Anti-CD54/Icam1 antibody [YN1/1.7.4] {Biotin} (STJA0001306)

STJA0001306

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Intercellular adhesion molecule 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	YN1/1.7.4
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	15894
Gene Symbol	Icam1
Uniprot ID	ICAM1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	