

Anti-Cd22 antibody [Cy34.1] (STJA0001313)

STJA0001313

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B-cell receptor CD22 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Cy34.1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12483
Gene Symbol	Cd22
Uniprot ID	CD22_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	