

Anti-Cd22 antibody [Cy34.1] {AF647} (STJA0001329)

STJA0001329

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B-cell receptor CD22 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Cy34.1
Concentration	0.5 mg/mL
Conjugation	AF647
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12483
Gene Symbol	Cd22
Uniprot ID	CD22_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	