

Anti-Cd209b antibody [22D1] (STJA0001330)

STJA0001330

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-CD209 antigen-like protein B is suitable for use in Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	FC, IHC-F, WB
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	22D1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	69165
Gene Symbol	Cd209b
Uniprot ID	C209B_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	