

Anti-CD127/IL7r antibody [A7R34] {PE/Cyanine7} (STJA0001338)

STJA0001338

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-7 receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A7R34
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16197
Gene Symbol	IL7r
Uniprot ID	IL7RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	