

Anti-CD272/Btla antibody [PK18.6] (STJA0001341)

STJA0001341

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-and T-lymphocyte attenuator is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PK18.6
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	208154
Gene Symbol	Btla
Uniprot ID	BTLA_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	