

Anti-CD272/Btla antibody [PK18.6] {Biotin} (STJA0001342)

STJA0001342

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-and T-lymphocyte attenuator is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PK18.6
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	208154
Gene Symbol	Btla
Uniprot ID	BTLA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	