

Anti-Cd28 antibody [37.51] (STJA0001343)

STJA0001343

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Hamster monoclonal antibody anti-T-cell-specific surface glycoprotein CD28 is suitable for use in Universal, Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Activ, Block, FC, IHC, IP
Host/Source	Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	37.51
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12487
Gene Symbol	Cd28
Uniprot ID	CD28_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	