

Anti-Cd28 antibody [37.51] (STJA0001344)

STJA0001344

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Hamster monoclonal antibody anti-T-cell-specific surface glycoprotein CD28 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications. |
| Applications | FC, IHC, IP |
| Host/Source | Hamster |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | 37.51 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG |
| Storage Instruction | Recommend storing at 2-8°C. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 12487 |
| Gene Symbol | Cd28 |
| Uniprot ID | CD28_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |