

Anti-Cd28 antibody [37.51] {FITC} (STJA0001349)

STJA0001349

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Hamster monoclonal antibody anti-T-cell-specific surface glycoprotein CD28 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	37.51
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12487
Gene Symbol	Cd28
Uniprot ID	CD28_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	