

Anti-Cd40 antibody [FGK4.5/FGK45] (STJA0001352)

STJA0001352

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FGK4.5/FGK45
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	21939
Gene Symbol	Cd40
Uniprot ID	TNR5_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	