

Anti-Cd40 antibody [FGK4.5/FGK45] {PE/Cyanine7} (STJA0001360)

STJA0001360

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FGK4.5/FGK45
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	21939
Gene Symbol	Cd40
Uniprot ID	TNR5_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	