

## Anti-CD122/IL2rb antibody [5H4] (STJA0001363)

STJA0001363

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Interleukin-2 receptor subunit beta is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
<b>Applications</b>	FC, IP
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5H4
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	Recommend storing at 2-8°C. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>	16185
<b>Gene Symbol</b>	IL2rb
<b>Uniprot ID</b>	IL2RB_MOUSE
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	